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Instructions For Use

Daily Disposable Contact Lenses

This instruction for use is intended for users of the contact lenses. Please read carefully and always follow your eye care professional's directions and the instructions.

EN

The information and instructions for safety and performance apply only to **Daily Wear Soft Contact Lenses**. For your eye health, it is important to wear your lenses as prescribed by your eye care professional. It is also important to keep your eye care professional fully aware of your medical history. Your eye care professional will tailor a total program of care based on your specific needs.

He or she will review with you all instructions for how to use the lens properly. The safe and effectiveness of contact lenses depends on proper use.

The light blue color will not make it more visible. The Soft Contact Lens with UV blocking may block UV light slightly, but cannot protect eyes as the UV protection equipment.

WARNING: UV absorbing contact lenses are NOT substitutes for protective UV absorbing eyewear such as UV absorbing goggles or sunglasses because they do not completely cover the eye and surrounding area. You should continue to use UV absorbing eyewear as directed.

Discard the lens if it becomes damaged or you reach the prescribed wearing time recommended by your eye care professional.

The Soft Contact Lens is daily wear contact lens prescribed for frequent replacement. They are different from the rigid contact as they are more soft and flexible.

WEARING RESTRICTIONS AND INDICATIONS :

The Spherical Contact Lenses Contact Lens is indicated as daily wear and single-use soft contact lens for the correction of refractive ametropia (myopia and hyperopia) in phakic or aphakic persons with non-diseased eyes who exhibit refractive astigmatism of 1.00D or less where the astigmatism does not interfere with visual acuity.

The Toric for Astigmatism Contact Lens is indicated as daily wear and single-use for the correction of ametropia (myopia and hyperopia) with astigmatism in aphakic and non-aphakic persons with non-diseased eyes and whose powers are from -12.00 to +8.00D and astigmatic corrections are from -0.75 to -2.25D.

The Multifocal for Presbyopia Contact Lens is indicated as daily wear and single-use for the correction of ametropia (myopia and hyperopia) and emmetropia with presbyopia in aphakic and non-aphakic persons with non-diseased eyes and whose powers are from -12.00 to +8.00D with add powers from +0.25 to +2.75D. The lenses may be worn by persons who exhibit astigmatism of 1.00D or less where the astigmatism does not interfere with visual acuity.

CONTRAINDICATIONS (REASONS NOT TO USE)

DO NOT USE Daily Wear Soft Contact Lenses when any of the following conditions exists:

- Acute or sub-acute inflammation or infection of the anterior chamber of the eye
- Any eye disease, injury or abnormality that affect the cornea, conjunctiva, or eyelids
- Severe insufficiency of lacrimal secretion (dry eye)
- Corneal hypoesthesia (reduced corneal sensitivity), if not-aphakic
- Any systemic disease that may affect the eye or be exacerbated by wearing contact lens
- Allergic reactions of ocular surface or adnexa that may be induced or exaggerated by the wearing of contact lenses or use of contact lens solutions
- Allergy to any ingredient, such as mercury or Thimerosal, in a solution which is be used to care for the **Daily Wear Soft Contact Lens**
- Any active corneal infection (bacterial, fungal, or viral)
- If eye become red or irritated

WARNINGS

INFECTION WITH POSSIBLE PERMANENT DAMAGE TO VISION COULD RESULT FROM THE FAILURE TO STRICTLY FOLLOW RECOMMENDED DIRECTIONS FOR USE.

PROBLEMS WITH CONTACT LENSES COULD RESULT IN SERIOUS INJURY TO THE EYE.

IT IS ESSENTIAL THAT YOU FOLLOW THE DIRECTIONS OF THE EYE CARE PROFESSIONAL AND ALL LABELING INSTRUCTIONS FOR PROPER USE OF CONTACT LENSES.

EYE PROBLEMS, INCLUDING CORNEAL ULCERS, CAN DEVELOP RAPIDLY AND LEAD TO LOSS OF VISION; IF YOU EXPERIENCE:

- Eye Discomfort
- Excessive Tearing
- Vision Changes
- Loss of Vision
- Eye Redness
- Or Other Eye Problems

YOU SHOULD BE INSTRUCTED TO IMMEDIATELY REMOVE THE LENSES, AND PROMPTLY CONTACT YOUR EYE CARE PROFESSIONAL.

- Daily wear lenses are not indicated for overnight wear, and you should be instructed not to wear lenses while sleeping.
- Please do not wear the contact lens continuously more than 12 hours .
- Always dispose of lenses when they are removed and have replacement lenses or glasses available.
- Consult the eye care professional regarding the use of contact lens in certain atmospheric or environmental conditions that can cause irritation to the eye.
- Consult the eye care professional or doctor if you are pregnant. According to personal physical condition, the dimension of cornea might have some change during pregnancy.
- Please note there are additional risks if you reuse of contact lens:
 - The chance of infection may increase and could lead to severe infection, vision loss or blindness.
 - You may experience blurred vision due to inappropriate cleaning of lens surface.
 - Lenses may become physically compromised by repeatedly removing and re-inserting lenses.
- Do not reuse old solution left in your lens case since solution reuse reduces effective lens disinfection and could lead to severe infection, vision loss or blindness.
- The lens is immersed in the packaging solution contents boric acid, to prevent microbiological contamination. If you might be pregnant, please talk to your doctor before wearing the contact lens as it contains boric acid.
- Borates, which include boric acid and sodium borate, are classified as Category 1B CMR substances when their concentration exceeds 0.1% by weight. However, these compounds are deemed safe for application when used in compliance with the provided label directions.
- Please flush it with clean saline solution before wearing the contact lens.

PRECAUTIONS

Eye care professionals should carefully instruct you about the following care regimen and safety precautions:

- If the lens sticks (stops moving) on the eye, you should be instructed to **immediately** consult your eye care professional.
- Always wash and rinse hands before handling lenses. Avoid contaminating lenses with make-up, aerosol sprays or wearing them near noxious fumes. It is best to put on lenses before putting on makeup.
- Avoid to scratches the lenses, causing distorted vision and/or injury of the eye, especially with fingernails.
- Carefully follow the handling, insertion, removal, and wearing instructions in you Instructions for the **Daily Wear Soft Contact Lens** and those prescribed by the eye care professional.
- Never wear lenses beyond the period recommended by the eye care professional.
- The lens has been sterilized and sealed in packaging until open for use. Do not attempt to re-sterilize the lens in any condition, even if it is unopened or unused. Inappropriate sterilization or re-sterilization may cause package to leak or damage the lens, which may lead to cornea scratches, irritation, sense of foreign object, infection or other possible abnormal conditions when wearing the damaged lens.
- Avoid all harmful or irritating vapors and fumes while wearing lenses.
- Ask the eye care professional about wearing lenses during sporting activities.
- Never use tweezers or other tools to remove lenses from the lens container unless specifically indicated for that use. Pour the lens into the hand.
- Always contact the eye care professional before using any medicine in the eyes.
- Always inform the employer of being a contact lens wearer. Some jobs may require use of eye protection equipment or may require that you not wear contact lenses.
- As with any contact lens, follow-up visits are necessary to assure the continuing health of your eyes. You should be instructed as to a recommend follow-up schedule.
- Never to share the lens with other users.

ADVERSE REACTIONS (PROBLEM AND WHAT TO DO)

- You should be informed that the following problems may occur:
- Eyes stinging, burning, itching (irritation), or other eye pain
 - Comfort is less than when lens was first placed on eye
 - Feeling that something is in the eye such as a foreign body or scratched area
 - Excessive watering (tearing) of the eyes
 - Unusual eye secretions
 - Redness of the eyes
 - Reduced sharpness of vision (poor visual acuity)
 - Blurred vision, rainbows, halos around objects
 - Sensitivity to light (photophobia)
 - Dry eyes


If you notice any of above, you should be instructed to:

- **Immediately remove lenses.**
- If the discomfort or problem stops, then look closely at the lens. If the lens is in any way damaged, Discard damaged lens. If you notice any problems, you should immediately remove your lens. If the problem or discomfort stops, discard the lens and place a new fresh lens on the eye.
- If your lens has dirt or foreign body on it, or the problem stops and the lens appears undamaged, you should thoroughly clean, rinse, and disinfect the lenses; then reinsert them. After reinsertion, if the problem continues, you should immediately remove the lenses and consult the eye care professional.

When any of the above problems occur, a serious condition such as infection, corneal ulcer, neovascularization, or iritis may be present. You should be instructed to **keep lens off the eye and seek immediate** professional identification of the problem and prompt treatment to avoid serious eye damage.

WEARING INSTRUCTION

The lens packages are individual and maintained sterilely (Sterilized using steam). DO NOT use if the sterile blister package is opened or damaged. If the lens appears damaged or deformed, DO NOT use it and use the next lens in the multi-pack.

- 1. Placing the Lens on the Eye**
- Wash and rinse your hands thoroughly and dry them thoroughly with a fiber-free towel, clean and dry before touching the lenses.
 - Verify that the lens is not turned inside out by placing it on your forefinger and checking its profile. The lens should appear a natural, curved, bowl-like shape.
 - Place the lens on the tip of the clean and dry right or left index finger. Place the middle finger of the same hand close to the lashes and pull the lower eyelid down.
 - Use the fingers of the other hand to lift the upper eyelid.
 - Place the lens directly on the eye (cornea) and gently spread your finger away from the lens.
 - Look down and slowly release your lower eyelid.
 - Look straight ahead and slowly release your upper eyelid.
 - Blink gently. If any problems discover, remove your lenses immediately and contact your eye care professional.

- 2. Centering the Lens**
- If the lens is not centered on the eye, to center a lens as follow one of the procedures below.
- a. Close your eyelids and gently massage the lens into place under the closed lids.
 - b. Use finger pressure on the edge of the upper lid or lower lid and gently manipulate the off-centered lens onto the cornea while the eye is open.

- 3. Removing the Lens**
- Wash, rinse, and dry your hands thoroughly.
- Caution:** Always be sure the lens is on the cornea before attempting to remove it. Simply follow the procedures recommended by the eye care professional.
- Pinch Method
 - a. Look up; slide the lens to the lower part of the eye using the forefinger.
 - b. Pinch the lens between the thumb and forefinger gently.
 - c. Remove the lens.
 - Remove the other lens by following the same procedure.
 - The proper method for disposing of used contact lenses is to throw the lenses in the trash bin.

REPORTING OF ADVERSE REACTIONS (Side Effects)

Any incidents should be reported to the manufacturer and/or its authorized representative and/or to your national authority.

It is essential that you follow your eye care professional's directions regarding this important step of building up your wear time.

SYMBOLS KEY

The following symbols may appear on the label or packaging.

Symbol	Definition	Symbol	Definition
	Consult Instructions for Use		Temperature limit
	Caution		Do not use if package is damaged
	Manufacturer		UV blocking
	Use-By date (expiration date)		CE-mark and Identification Number of Notified Body
	Batch code		Authorized representative in the European Community
	Sterilized using steam		Medical device
	Diameter (the unit is millimeter)		Unique device identifier
	Base curve (the unit is millimeter)		Contains hazardous substances
	Diopter (lens power)		Do not re-use
	Cylinder		Date of manufacture
	Axis		Single sterile barrier system
	Highest near addition that can be corrected		Importer
	"low" near ADD lens (up to +1.25D)		Do not resterilize
	"medium" near ADD lens (up to +2.00D)		"high" near ADD lens (above +2.75D)

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Emergo Europe B.V.